SARS-CoV-2 Ag Rapid Test

Package Insert

REF VCD05-01-011/ VCD05-01-012/ VCD05-01-013

English

PRINCIPLE AND INTENDED USE

VivaDiag™ SARS-CoV-2 Ag Rapid Test is for the rapid, qualitative detection of the nucleocapsid protein antigen from SARS-CoV-2 in human nasal swab, oropharyngeal swab or nasopharyngeal swab specimen. The test is for *in vitro* diagnostic use only. For professional use only. It is intended for clinical laboratories and healthcare professional use only as point-of-care testing. Not for at-home testing.

VivaDiag™ SARS-CoV-2 Ag Rapid Test is based on immunoassay technology. Each test device has one line of anti-SARS-CoV-2 monoclonal antibody on the detection line (T line) and one line of anti-mouse IgG polyclonal antibody on the quality control line (C line). When extracted specimen is added to the specimen well, it will react with the labeled antibody to form a complex, the mixture then migrates through the membrane by capillary action and interacts with the coated anti-SARS-CoV-2 monoclonal antibody on the detection line. If the specimen contains SARS-CoV-2 antigen, the detection line will appear purplish-red indicating the SARS-CoV-2 antigen is positive. Otherwise, the test result will be negative. The test device also contains a quality control line C which should appear purplish-red for all valid tests. If the quality control line C does not appear, the test result will be invalid even if the detection line appears.

COMPOSITION

Each test kit contains test devices, sealed pouches (prefilled with 300 μ L extraction solution), extraction tubes, extraction tube tips, tube stand, sterile swabs and package insert.

Materials required but not provided: timer.

STORAGE AND HANDLING

• Store the test kit in a cool, dry place between 2-30°C. Keep away from light. Exposure to temperature and / or humidity outside the specified conditions may cause inaccurate results.

• Do not freeze or refrigerate. Use the test kit at temperatures between 15-30°C.

• Use the test kit between 10-90% humidity.

• Do not use the test kit beyond the expiration date (printed on the foil pouch and box).

Note: All expiration dates are printed in Year-Month-Day format. 2022-06-18 indicates June 18, 2022.

WARNINGS, PRECAUTIONS AND LIMITATIONS

• Results from SARS-CoV-2 antigen testing should not be used as the sole basis to diagnose or exclude SARS-CoV-2 infection or to inform infection status.

• Negative results do not rule out SARS-CoV-2 infection, particularly in those who have been in contact with the virus. Follow-up testing with a molecular diagnostic and / or CT should be considered to rule out infection in these individuals.

• Positive results may be due to present infection with SARS-coronavirus strains, see "cross-reactivity" for details. Follow-up testing with a molecular diagnostic and / or CT should be considered to confirm the testing result.

• For *in vitro* diagnostic use only.

• Not for at-home testing.

• Further molecular diagnostic and / or CT is recommended to identify the actual physical situation.

• Do not open the foil pouch of the test device exposing it to the ambient environment until the test device is ready for immediate use.

• Do not use any damaged test device or material.

• Do not reuse the test device.

• Handle extraction solution with caution, do not contact with eyes or skin. If spilled on eyes or skin, wash thoroughly with water.

• Do not test use kit beyond the expiration date.


• Specific training or guidance is recommended if operators are not experienced with specimen collection and handling procedures.

• Only use nasal swab, oropharyngeal swab or nasopharyngeal swab as specimen. Follow the package insert to obtain accurate results.

• Wear protective gears such as laboratory coats, disposable gloves, and eye protection when specimens are collected and evaluated.

• Wash hands thoroughly after handling.

• All parts of kit are considered biohazardous and can potentially transmit infectious diseases from blood-borne pathogens, even after you have performed cleaning and disinfection. Follow proper precautions and all local regulations when disposing of the used test kits.

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Lowest Concentration with Uniform Positivity per Analyte	6.75x10 ⁵ TCID ₅₀ /mL
Limit of detection (LoD) per inactivated Virus Culture	6.75x10 ⁵ TCID ₅₀ /mL

2. Clinical Sensitivity/Clinical Specificity

A total of 552 specimens were tested using the VivaDiag™ SARS-CoV-2 Ag Rapid Test. These specimens were obtained by nasal swabs from symptomatic patients. The performance of the VivaDiag™ SARS-CoV-2 Ag Rapid Test was compared to a commercialized molecular assay.

Table Summary of sensitivity/specificity of the VivaDiag™ SARS-CoV-2 Ag Rapid Test compared to PCR.

VivaDiag™ SARS-CoV-2 Ag Rapid Test	PCR		
	Positive	Negative	Total
Positive	115	0	115
Negative	6	431	437
Total	121	431	552
Sensitivity	95.04%		
(115/121, 95%CI, 89.60%~97.71%)			
Specificity	100%		
(431/431, 95%CI, 99.12%~100%)			
Accuracy	98.91%		
(546/552, 95%CI, 97.65%~99.50%)			

The VivaDiag™ SARS-CoV-2 Ag Rapid Test showed a clinical sensitivity of 95.04%.

The VivaDiag™ SARS-CoV-2 Ag Rapid Test showed a clinical specificity of 100%.

The VivaDiag™ SARS-CoV-2 Ag Rapid Test showed a clinical accuracy of 98.91%.

CROSS-REACTIVITY

1. Cross-Reactivity: there was no cross-reaction with potential cross-reactive substances except SARS-coronavirus.

1) cross-reaction with SARS-coronavirus.

Virus	Strain	Concentration
SARS-coronavirus	Urbani	1X10 ⁶ PFU/mL

2) no cross-reaction with potential cross-reactive substances

Virus/Bacteria/Parasite	Strain	Concentration Range
	H1N1	
Influenza A	H3N2	
	H5N1	
	H7N9	
Influenza B	NA	
	Type1	
	Type2	
Adenovirus	Type3	
	Type5	
	Type7	
Respiratory syncytial virus	Type55	1X10 ⁴ ~1X10 ⁶ TCID ₅₀ /mL
	Type A	
	Type B	
Coronavirus	229E	
	OC43	
	NL63	
MERS-Coronavirus	Florida/USA-2_Saudi Arabia,2014	
Parainfluenza virus	Type1	
	Type2	
	Type3	
	Type4	

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