

HIV Real-TM Qual

Handbook

Real Time PCR kit for qualitative detection of
Human Immunodeficiency Virus (HIV)

REF R-V0-100FRT

REF TR-V0-100FRT

▽ 100

NAME

HIV Real-TM Qual

INTENDED USE

kit **HIV Real-TM Qual** is a Real-Time test for the qualitative detection of Human Immunodeficiency Virus (HIV) and simultaneous detection of a HIV-specific Internal Control (IC), by dual color detection.

PRINCIPLE OF ASSAY

kit **HIV Real-TM Qual** is a Real-Time test for the detection of Human Immunodeficiency Virus in human plasma. HIV RNA is extracted from plasma, amplified using real time amplification and detected using fluorescent reporter dye probes specific for HIV or HIV IC.

Internal Control (IC) serves as an amplification control for each individually processed specimen and to identify possible inhibition. IC is detected in a channel other than the HIV RNA. Monitoring the fluorescence intensities during Real Time allows the detection of the accumulating product without having to re-open the reaction tube after the real time amplification.

MATERIALS PROVIDED

Module No.1: Real Time PCR kit (R-V0-100FRT)

Part N° 2 – “Controls”

- **HIV Rec RNA Pos Control***, 4 x 0,01 ml;
- **Negative Control****, 4 x 0,5 ml;
- **HIV Rec IC (Internal Control)*****, 4 x 0,13 ml.

Part N° 3 – “HIV Real-TM Qual”:

- **DTT**, 4 tubes.
- **RT-PCR-mix-1-TM HIV**, 4 x 0,3 mL.
- **RT-PCR-mix-2-TM**, 4 x 0,2 mL.
- **Hot Start TaqF Polymerase**, 4 x 0,02 mL
- **M-MLV Revertase**, 4 x 0,01 mL;
- **TE-buffer**, 4 x 0,07 mL
- **Standard HIV**
 - **QS3 HIV**, 4 x 0,025 mL;
- **Standard IC**
 - **QS3 IC**, 4 x 0,025 mL;

The kit allows to detect 100 samples in 25 µl reaction mix

* *must be used in the isolation procedure as Positive Control of Extraction (when using SaMag Viral Nucleic Acids extraction kit always add 10 µl of Positive Control of Extraction regardless of the extraction starting volume. For example add 10 µl of Pos Control and 390 µl of negative control if the extraction starting volume is 400 µl).*

** *must be used in the isolation procedure as Negative Control of Extraction.*

*** *add 5 µl of Internal Control during the RNA isolation directly to the sample/lysis mixture*

Module No.2: Real Time PCR kit (TR-V0-100FRT)

Part N° 1 – “Ribo-Sorb-100”:

- **Lysis Solution**, 2 x 22,5 ml;
- **Washing Solution**, 2 x 20,0 ml;
- **Sorbent**, 2 x 1,25 ml;
- **RNA-eluent**, 10 x 0,5 ml

Part N° 2 – “Controls”

- **HIV Rec RNA Pos Control***, 4 x 0,01 ml;
- **Negative Control****, 4 x 0,5 ml;
- **HIV Rec IC (Internal Control)*****, 4 x 0,13 ml.

Part N° 3 – “HIV Real-TM Qual”:

- **DTT**, 4 tubes.
- **RT-PCR-mix-1-TM HIV**, 4 x 0,3 mL.
- **RT-PCR-mix-2-TM**, 4 x 0,2 mL.
- **Hot Start TaqF Polymerase**, 4 x 0,02 mL
- **M-MLV Revertase**, 4 x 0,01 mL;
- **TE-buffer**, 4 x 0,07 mL
- **Standard HIV**
 - **QS3 HIV**, 4 x 0,025 mL;
- **Standard IC**
 - **QS3 IC**, 4 x 0,025 mL;

The kit allows to detect 100 samples in 25 µl reaction mix

* *must be used in the isolation procedure as Positive Control of Extraction.*

** *must be used in the isolation procedure as Negative Control of Extraction.*

*** *add 5 µl of Internal Control during the RNA isolation directly to the sample/lysis mixture*

MATERIALS REQUIRED BUT NOT PROVIDED

Zone 1: sample preparation:

- Biological cabinet
- Desktop microcentrifuge for “eppendorf” type tubes (RCF max. 16,000 x g); Eppendorf 5415D or equivalent
- 60°C ± 2°C dry heat block
- Vortex mixer
- Pipettors (capacity 5-40 µl; 40-200 µl; 200-1000 µl) with aerosol barrier
- 1,5 ml polypropylene sterile tubes (Sarstedt, QSP, Eppendorf)
- Disposable gloves, powderless
- Biohazard waste container
- Refrigerator, Freezer
- 70% Ethanol (freshly prepared mixture of reagent grade 96% ethanol and distilled water)
- Acetone

Zone 2: RT and amplification:

- Real Time Thermal cycler
- Tubes (0,2 ml)
- Workstation
- Pipettors (capacity 0,5-10 µl; 5-40 µl) with aerosol barrier
- Tube racks

STORAGE INSTRUCTIONS

Part N° 1 – “**Ribo-Sorb-100**” must be stored at 2-8°C

Part N° 2 – “**Controls**” must be stored at 2-8°C.

Part N° 3 – “**HIV Real-TM Qual**” must be stored at -20°C.

The kits can be shipped at 2-8°C for 3-4 days but should be stored at 2-8°C and -20°C immediately on receipt.

STABILITY


HIV Real-TM Qual Test is stable up to the expiration date indicated on the kit label. The product will maintain performance through the control date printed on the label. Exposure to light, heat or humidity may affect the shelf life of some of the kit components and should be avoided. Repeated thawing and freezing of these reagents should be avoided, as this may reduce the sensitivity. Components stored under conditions other than those stated on the labels may not perform properly and may adversely affect the assay results.

QUALITY CONTROL

In accordance with Sacace's ISO 13485-Certified Quality Management System, each lot is tested against predetermined specifications to ensure consistent product quality.

WARNINGS AND PRECAUTIONS

The user should always pay attention to the following:

-  Lysis Solution contains guanidine thiocyanate*. Guanidine thiocyanate is harmful if inhaled, or comes into contact with skin or if swallowed. Contact with acid releases toxic gas. (Xn; R: 20/21/22-36/37/38; S: 36/37/39).
- Use sterile pipette tips with aerosol barriers and use new tip for every procedure.
- Store extracted positive material (samples, controls and amplicons) away from all other reagents and add it to the reaction mix in a separate area.
- Thaw all components thoroughly at room temperature before starting an assay.
- When thawed, mix the components and centrifuge briefly.
- Use disposable gloves, laboratory coats and eye protection when handling specimens and reagents. Thoroughly wash hands afterwards.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in laboratory work areas.
- Do not use a kit after its expiration date.
- Dispose of all specimens and unused reagents in accordance with local authorities' regulations.
- Specimens should be considered potentially infectious and handled in a biological cabinet in accordance with appropriate biosafety practices.
- Clean and disinfect all sample or reagent spills using a disinfectant such as 0.5% sodium hypochlorite, or other suitable disinfectant.
- Avoid sample or reagent contact with the skin, eyes, and mucous membranes. If skin, eyes, or mucous membranes come into contact, rinse immediately with water and seek medical advice immediately.
- Material Safety Data Sheets (MSDS) are available on request.
- Use of this product should be limited to personnel trained in the techniques of DNA amplification.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.



Some components of this kit contain sodium azide as a preservative. Do not use metal tubing for reagent transfer.

PRODUCT USE LIMITATIONS

Use of this product should be limited to personnel trained in the techniques of DNA amplification (UNI EN ISO 18113-2:2012). Strict compliance with the user manual is required for optimal PCR results. Attention should be paid to expiration dates printed on the box and labels of all components. Do not use a kit after its expiration date.

SAMPLE COLLECTION, STORAGE AND TRANSPORT

HIV Real-TM Qual can analyze RNA extracted from:

- *plasma* collected blood in ACD or EDTA tubes;

Specimens can be stored at +2-8°C for no longer than 12 hours, or freeze at -20°C to -80°C.

Transportation of clinical specimens must comply with country, federal, state and local regulations for the transport of etiologic agents.

RNA ISOLATION

The following kits are recommended:

- ⇒ **Ribo-Sorb** (Sacace, REF K-2-1)
- ⇒ **Ribo-Virus** (Sacace, REF K-2-C)
- ⇒ **SaMag Viral Nucleic Acids Extraction kit** (Sacace, REF SM003)

Please carry out RNA extraction according to the manufacture's instruction. Add 5 µl of Internal Control during RNA isolation procedure directly to the sample/lysis mixture.

INTERNAL CONTROL (HIV Rec IC)

HIV Rec IC is a qualitative Internal Control and represents recombinant RNA-containing-structure which carried through all steps of analysis from nucleic acid extraction to PCR amplification-detection. The presence of HIV Rec IC allows to monitor the extraction procedure and to check possible PCR inhibition.

SPECIMEN AND REAGENT PREPARATION (Only for **Kit TR-V0-100FRT**)

To increase the sensitivity of test it is recommended to ultracentrifuge plasma samples (1 ml) at 24,000 x g at 4°C for 60 minutes. Discard the supernatant (900 µl) and use the remaining 100 µl pellet for RNA extraction.

1. **Lysis Solution** and **Washing Solution** should be warmed up to 60–65°C until disappearance of ice crystals.
2. Prepare 70% Ethanol.
3. Prepare the required quantity of 1,5 ml polypropylene tubes including one tube for **Negative Control of Extraction** and two tubes for **Positive Controls of Extraction**.
4. Add to each tube **5 µl of Internal Control** and **450 µl Lysis Solution**.
5. Vortex thawed patient plasma specimens for 5 sec.
6. Add **450 µl of Lysis Solution with IC** to each of labeled tubes.
7. Add **100 µl of Samples** to the appropriate tube.
8. Prepare Controls as follows:
 - add **100 µl of C– (Negative Control)** to labeled *Cneg*.
 - add **90 µl of C– (Negative Control)** and **10 µl of HIV Rec Pos** to the tube labeled *Cpos1*.
9. Vortex the tubes and centrifuge for 7-10 sec.
10. Vortex vigorously **Sorbent** and add **25 µl** to each tube.
11. Vortex for 5-7 sec and incubate all tubes for 10 min at room temperature. Vortex periodically.
12. Centrifuge all tubes for 1 min at 10000g and using a micropipette or vacuum aspirator with a plugged aerosol barrier tip, carefully remove and discard supernatant from each tube without disturbing the pellet. Change tips between tubes.
13. Add **400 µl of Washing Solution** to each tube. Vortex vigorously until the sorbent is completely resuspended, centrifuge for 1 min at 10000g and using a micropipette or vacuum aspirator with a plugged aerosol barrier tip, carefully remove and discard supernatant from each tube without disturbing the pellet. Change tips between tubes.
14. Add **500 µl of Ethanol at 70%** to each tube. Vortex vigorously until the sorbent is completely resuspended, centrifuge for 1 min at 10000g and using a micropipette or vacuum aspirator with a plugged aerosol barrier tip, carefully remove and discard supernatant from each tube without disturbing the pellet. Change tips between the tubes.
15. Repeat step 14.
16. Add **500 µl of Acetone** to each tube. Vortex vigorously and centrifuge for 1 min at 10000g and using a micropipette or vacuum aspirator with a plugged aerosol barrier tip, carefully remove and discard supernatant from each tube without disturbing the pellet. Change tips between the tubes.
17. Incubate all tubes with open cap for 10 min at 56°C.

18. Resuspend the pellet in **50 µl** of **RNA-eluent**. Incubate for 10 min at 56°C and vortex periodically.

Centrifuge the tubes for 2 min at maximum speed (12000-16000 g). The supernatant contains RNA/DNA ready for amplification. The Real Time amplification can be performed on the same day of extraction.

RT AND AMPLIFICATION

1. Thaw one set of reagents, vortex and centrifuge briefly the tubes.
2. Prepare 0,2 µl tubes.
3. Prepare **Reaction Mix**: add into the tube with **DTT 300 µl** of **RT-PCR-mix-1 HIV**, **200 µl** of **RT-PCR-mix-2**, **20 µl** of **Hot Start TaqF Polymerase** and **7,5 µl** of **M-MLV Revertase**. Vortex thoroughly and centrifuge briefly.

*(If it is necessary to test less than 25 samples add into the tube with **DTT**, **300 µl** of **RT-PCR-mix-1**, **200 µl** of **RT-PCR-mix-2** and vortex for at least 5-10 seconds. This mix is stable for 1 month at -20°C. Add for each sample (N) in the new sterile tube **12,5*N µl** of mix, **0,5 *N µl** of **TaqF Polymerase** and **0,25 *N µl** of **M-MLV**.)*

4. Add **12,5 µl** of **Reaction Mix** into each tube.
5. Add **12,5 µl** of **extracted RNA** sample to the appropriate tube with Reaction Mix and mix by pipetting

If the Ribo-Sorb isolation kit is used as a RNA extraction kit, re-centrifuge all the tubes with extracted RNA for 2 min at maximum speed (12000-16000 g) and take carefully supernatant.

N.B. don't disturb the pellet, sorbent inhibit reaction!

6. Prepare for each run 2 standards and 1 negative control:
 7. add **12,5 µl** of **QS3 HIV** into HIV Pos labeled tube;
 8. add **12,5 µl** of **QS3 IC** into IC Pos labeled tube;
 9. add **12,5 µl** of **TE-buffer** to the tube labeled Negative Control;

Amplification

1. Create a temperature profile on your instrument as follows:

Step	Rotor-type Instruments ¹			Plate- or modular type Instruments ²		
	Temperature, °C	Time	Repeats	Temperature, °C	Time	Repeats
1	50	30 min	1	50	30 min	1
2	95	15 min	1	95	15 min	
2	95	20 s	5	95	20 s	5
	52	30 s		52	30 s	
	72	30 s		72	30 s	
3	95	20 s	40	95	20 s	42
	55	30 s <i>fluorescent signal detection</i>		55	40 s <i>fluorescent signal detection</i>	
	72	30 s		72	30 s	

¹ For example Rotor-Gene™ 3000/6000/Q (Corbett Research, Qiagen)

² For example, SaCycler-96™ (Sacace), CFX/iQ5™ (BioRad); Mx3005P™ (Agilent), ABI® 7300/7500/StepOne Real Time PCR (Applied Biosystems)

Fluorescence is detected at the 2nd step of Cycling 2 stage (55 °C) in FAM/Green and JOE/Yellow/Hex/Cy3 fluorescence channels.

HIV RNA is detected on the JOE(Yellow)/HEX/Cy3, IC RNA on the FAM (Green) channel

INSTRUMENT SETTINGS

Rotor-type instruments

Channel	Calibrate/Gain Optimisation...	Threshold	More Settings/ Outlier Removal	Slope Correct
FAM/Green	<i>from 5 FI to 10 FI</i>	0.03	10 %	ON
JOE/Yellow	<i>from 5 FI to 10 FI</i>	0.03	5 %	ON

Plate-type instruments

The threshold line should cross only sigmoid curves of signal accumulation of positive samples and should not cross the baseline; otherwise, the threshold level should be raised. Set the threshold at a level where fluorescence curves are linear and do not cross curves of the negative samples.

DATA ANALYSIS

The fluorescent signal intensity is detected in two channels:

The signal from the Internal Control amplification product is detected in the FAM fluorescence channel;

The signal from the *HIV* RNA amplification product is detected in the JOE channel.

Interpretation of results

The results are interpreted by the Instrument software by the crossing (or not-crossing) of the fluorescence curve with the threshold line and are shown as the presence (or absence) of Ct (threshold cycle) in the result grid.

The result of the analysis is considered reliable only if the results for Positive and Negative Controls of Amplification as well as Negative Control of Extraction are correct.

1. The sample is considered to be **positive** if its Ct value detected in the result grid in the channel for Positive Control is less than the specified boundary Ct value.
2. The sample is considered to be **negative** if its Ct value is not detected in the result grid in the channel for Positive Control (the fluorescence curve does not cross the threshold line) and if the Ct value determined in the results grid in the channel for IC does not exceed the specified boundary Ct value.

Boundary value of the cycle threshold, Ct

Sample	Ct boundary value Plate/Rotor Type Instrument	
	FAM/Green	JOE/Yellow/HEX
C-	38	-
PCE	38	40
NCA	-	-
QS3 HIV	-	40
QS3 IC	38	-
Clinical samples	38	40

QUALITY CONTROL PROCEDURE

A defined quantity of Internal Control (IC) is introduced into each sample and control at the beginning of sample preparation procedure in order to control the extraction process of each individual sample and to identify possible reaction inhibition.

A negative control of extraction (NCE), negative amplification control (NCA), positive amplification control (C+) are required for every run to verify that the specimen preparation, the amplification and the detection steps are performed correctly.

If the controls are out of their expected range (see table Results for Controls), all of the specimens and controls from that run must be processed beginning from the sample preparation step.

PERFORMANCE CHARACTERISTICS

Analytical specificity

The analytical specificity of the primers and probes was validated with negative samples. They did not generate any signal with the specific **HIV** primers and probes. The specificity of the kit **HIV Real-TM Qual** was 100%. The potential cross-reactivity of the kit **HIV Real-TM Qual** was tested against the group control. It was not observed any cross-reactivity with other pathogens.

Analytical sensitivity










The kit **HIV Real-TM Qual** allows to detect **HIV** RNA in 100% of the tests with a sensitivity of not less than 50 copies/ml . The detection was carried out on the control standard and its dilutions by negative plasma using the “Magno-Sorb” extraction kit (Sacace REF K-2-16/1000).

Target region: pol gene

TROUBLESHOOTING

1. Weak or no signal of the IC (FAM/Green channel) for the Negative Control of extraction.
 - The PCR was inhibited.
 - ⇒ Make sure that you use a recommended RNA extraction method and follow to the manufacturer's instructions.
 - ⇒ Re-centrifuge all the tubes before pipetting of the extracted RNA for 2 min at maximum speed (12000-16000 g) and take carefully supernatant. Don't disturb the pellet, sorbent inhibit reaction.
 - The reagents storage conditions didn't comply with the instructions.
 - ⇒ Check the storage conditions
 - Improper RNA extraction.
 - ⇒ Repeat analysis starting from the RNA extraction stage
 - The PCR conditions didn't comply with the instructions.
 - ⇒ Check the PCR conditions and select for the IC detection the fluorescence channel reported in the protocol.
 - The IC was not added to the sample during the pipetting of reagents.
 - ⇒ Make attention during the RNA extraction procedure.
2. Weak or no signal of the Positive Control.
 - The PCR conditions didn't comply with the instructions.
 - ⇒ Check the amplification protocol and select the fluorescence channel reported in the manual.
3. JOE(Yellow)/HEX/Cy3 signal with Negative Control of extraction.
 - Contamination during RNA extraction procedure. All samples results are invalid.
 - ⇒ Decontaminate all surfaces and instruments with sodium hypochlorite and ethanol.
 - ⇒ Use only filter tips during the extraction procedure. Change tips between tubes.
 - ⇒ Repeat the RNA extraction with the new set of reagents.
4. Any signal with Negative Control of PCR (DNA-buffer).
 - Contamination during PCR preparation procedure. All samples results are invalid.
 - ⇒ Decontaminate all surfaces and instruments with sodium hypochlorite and ethanol or special RNA decontamination reagents.
 - ⇒ Pipette the Positive control at last.
 - ⇒ Repeat the PCR preparation with the new set of reagents.

KEY TO SYMBOLS USED

	List Number		Caution!
	Lot Number		Contains sufficient for <n> tests
	Expiration Date		Version
	Store at	NCA	Negative Control of Amplification
	Manufacturer	NCE	Negative control of Extraction
	Consult instructions for use	C+	Positive Control of Amplification
		IC	Internal Control

- * SaCycler™ is a registered trademark of Sacace Biotechnologies
- * CFX™ and iQ5™ are registered trademarks of Bio-Rad Laboratories
- * Rotor-Gene™ is a registered trademark of Qiagen
- * MX3005P® is a registered trademark of Agilent Technologies
- * ABI® is a registered trademark of Applied Biosystems



Sacace Biotechnologies Srl
 via Scalabrini, 44 – 22100 – Como – Italy Tel +390314892927 Fax +390314892926
 mail: info@sacace.com web: www.sacace.com

