

CDV/CAV Ag

General Information

Purpose : Detection of Canine Distemper virus and Adeno virus antigen

Principle : Immunochromatographic assay


Specimen : CDV Ag - Conjunctiva, CAV Ag - Nasal fluid

Materials provided

- Rapid CDV/CAV test device
- Assay diluents tube
- Disposable swab
- Disposable dropper

Special Features

- Simultaneous detection of CDV and CAV with one specimen preparation
- One-step testing procedure : Quick & Accurate
- High sensitivity and specificity




Gentaur Group
Europe, UK & USA
info@gentaur.com



Order / Buy Now

Test Procedure



Ordering Information

Cat. No.	Description	Type	Packing size
173-RG11-08	CDV/CAV Ag	Device	10 Tests/Kit