


Zeptomatrix NATtrol GI Panels Compatibility with BioFire® FilmArray® GI Panel



info@gentaur.com

BIOFIRE® FILMARRAY®GI Panel	NATGIP-BIO	Notes	NATGIP-ABC	Notes
Campylobacter (jejuni, coli & upsaliensis)	✓	w/o upsaliensis	✓	w/o coli & upsaliensis
Clostridium difficile (Toxin A/B)	✓		✓	
Plesiomonas shigelloides	✓		X	
Salmonella	✓	enterica & typhimurium	X	
Yersinia enterocolitica	✓		✓	
Vibrio (parahaemolyticus, vulnificus, & cholerae)	✓	w/o parahaemolyticus & vulnificus	✓	w/o vulnificus
Vibrio cholerae	✓		✓	
E. coli O157	✓		X	
Enteroaggregative E. coli (EAEC)	✓		✓	
Enteropathogenic E. coli (EPEC)	✓		✓	
Enterotoxigenic E. coli (ETEC) lt/st	✓		✓	
Shiga-like toxin-producing E. coli (STEC) stx1/stx2 E. coli O157	X		X	
Shigella/Enteroinvasive E. coli (EIEC)	✓	S. sonnei	✓	S. sonnei
Adenovirus F 40/41	✓	type 41	✓	type 41
Astrovirus	✓		X	
Norovirus GI/GII	✓		✓	
Rotavirus A	✓	Strain Wa	✓	Strain Wa
Sapovirus (I,II, IV, and V)	✓		X	
Cryptosporidium	✓	C. parvum, strain Iowa	✓	C. parvum, strain Iowa
Cyclospora cayetanensis	✓		X	
Entamoeba histolytica	✓	strain DS4-868	✓	strain DS4-868
Giardia lamblia	✓	strain H3	X	



[NATGIP-BIO | NATtrol™ GI Panel \(22 x 0.25mL; 4 x 0.85 mL\)](#)
[NATGIP-ABC | NATtrol™ GI Panel \(17 x 0.25mL; 2 x 0.85 mL\)](#)

[Request Quotation](#)