


# Zeptomatrix NATtrols Compatible with BioFire® FilmArray® Respiratory 2.1 plus Panel



info@gentaur.com



BIOFIRE® Respiratory 2.1 plus Panel	NATRPC2-BIO	Notes	NATRVP2-QIA	Notes	NATRVP-IDI	Notes
Adenovirus	✓	types 1, 3, 31	✓	type 3	✓	type 3
Coronavirus 229E	✓		✓		✓	
Coronavirus HKU1	✓		✓		✓	
Coronavirus OC43	✓		✓		✓	
Coronavirus NL63	✓		✓		✓	
Middle East Respiratory Syndrome CoronaVirus (Mers-CoV)	X		X		X	
Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2)	X		✓	strain USA-WA1/2020	X	
Human Metapneumovirus	✓	strain Peru6-2003	✓	strain Peru6-2003	✓	strain Peru6-2003
Human Rhinovirus/Enterovirus	✓	Rhinovirus strain 1A	✓	Rhinovirus strain 1A	✓	Rhinovirus strain 1A
Influenza A	X	no general Influenza A	X	no general Influenza A	X	no general Influenza A
Influenza A/H1	✓	type H1N1 strain A/New Caledonia/20/99	✓	type H1N1 strain A/New Caledonia/20/99	✓	type H1N1 strain A/New Caledonia/20/99
Influenza A/H1-2009	✓	strain A/NY/02/09	✓	strain A/NY/02/09	✓	strain A/NY/02/09
Influenza A/H3	✓	strain A/Brisbane/10/07	✓	strain A/Brisbane/10/07	✓	strain A/Brisbane/10/07
Influenza B	✓	strain B/Florida/02/06	✓	strain B/Panama/45/90	✓	strain B/Florida/02/06
Parainfluenza 1	✓		✓		✓	
Parainfluenza 2	✓		✓		✓	
Parainfluenza 3	✓		✓		✓	
Parainfluenza 4	✓		✓		✓	
RSV	✓	type A, strain 2006 isolate	✓	type A	✓	type A
Bordetella pertussis	✓	strain A639	✓	strain A639	✓	strain A639
Bordetella parapertussis	✓	strain A747	X		X	
Chlamydomphila pneumoniae	✓	strain CWL-029	✓	strain CWL-029	✓	strain CWL-029
Mycoplasma pneumoniae	✓	strain M129	✓	strain M129	✓	strain M129

NATRPC2-BIO | Respiratory Panel 2 (RP2) Controls (Ea)  
 NATRVP2-QIA | Respiratory Verification Panel 2 (21 x 0.25mL) (EA)  
 NATRVP-IDI | Respiratory Verification Panel (RVP) (20 X 0.6 mL)

[Request Quotation](#)