

Swine

Real-Time PCR Kit

PRRSV, PCV2, CSFV detection Kit



PCR

Background

***PCV-2 Real -Time Detection Kit and *CSFV Real -Time Detection Kit from BIONOTE are** useful for detection the infection of Porcine Circovirus or Classical swine fever virus. The Kit can be exactly performed to detect DNA of Porcine Circovirus or Classical swine fever virus, so it can be used for both qualitative and quantitative analysis. It works most of Real-Time PCR apparatuses of block and capillary tubes

***PRRSV Real -Time Detection Kit from BIONOTE is** useful for **differential** detection the infection of EU/US strain of Porcine Reproductive and Respiratory Syndrome Virus. The Kit can be exactly performed to detect DNA of Porcine Reproductive and Respiratory Syndrome Virus (EU/US), so it can be used for both qualitative and quantitative analysis. It works most of Real-Time PCR apparatuses of block and capillary tubes


PCV-2 or CSFV Real-Time Detection Kit

Contents	Volume
Standard 1 (1x10 ⁴ copies/μl)	200μl x1
Standard 2 (1x10 ⁵ copies/μl)	200μl x1
Standard 3 (1x10 ⁶ copies/μl)	200μl x1
Standard 4 (1x10 ⁷ copies/μl)	200μl x1
Detection Solution PCV-2	500μl x1
Nuclease Free Water	1500μl x1
2X Master Mix	1000μl x1
Rox Reference Dye	40μl x1

PRRSV Real-Time Detection Kit

Contents	Volume
Standard 1 (1x10 ⁴ copies/μl)	200μl x1
Standard 2 (1x10 ⁵ copies/μl)	200μl x1
Standard 3 (1x10 ⁶ copies/μl)	200μl x1
Standard 4 (1x10 ⁷ copies/μl)	200μl x1
Detection Solution PRRS EU	460μl x1
Detection Solution PRRS US	460μl x1
Nuclease Free Water	1500μl x1
2X Enzyme buffer	1000μl x1
Enzyme Mix	40μl x1
Rox Reference Dye	40μl x1

Swine Real-Time Detection Kits Procedure



List of PCV and PRRSV Real-Time Detection Kit

Product	Diagnosis	Cat.No.
PCV-2 Real-Time Detection Kit	Porcine Circovirus	PD64-03
PRRSV Real-Time Detection Kit	Porcine Reproductive and Respiratory Syndrome Virus (Differential diagnosis of EU/US strain)	PD64-04
CSFV Real-Time Detection Kit	Classical swine fever virus	PD64-05